

R ROBOTIC DEVICES


AMADEO®

AMADEO® is the pioneer in hand and finger rehabilitation. AMADEO® helps to train various aspects of the **grasping movement**, offers targeted **finger extension** training, and enables to treat patients with **spasticity**.



DIEGO®

DIEGO® is a unique therapy device for shoulder and arm rehabilitation. The **intelligent gravity support** allows **active training** at the patient's individual capacity in a **three-dimensional space**.



OMEGO® PLUS

OMEGO® Plus is a **robotic therapy system** that trains multiple aspects needed for walking (motor function, range of motion, strength, endurance, balance, coordination, and cognition). The **highly adjustable** therapy chair makes OMEGO® Plus accessible for patients in all phases of rehabilitation.



LEXO®

LEXO® is an **end-effector robotic gait training** system, allowing free joint movement and encouraging the patient's self-activity. Easy handling and fast setup allow for **efficient use of therapy time**.



TyroS

The TyroS software is the heart of all Tyromotion devices, combining technology, know-how, and therapeutic programs. The software offers comprehensive patient view across all devices and encourages patients with its motivational therapeutic gaming environment.

For more information visit our website or contact us

SENSORY DEVICES

MYRO®

MYRO® is a sensor-based therapy device for motor rehabilitation of the upper extremities and **cognitive therapy**. MYRO® offers **force detection**, integration of **real objects** and a large variety of **gamification-based** therapy programs.



PABLO®

PABLO® is a sensor-based rehabilitation device for impaired motor function. PABLO® allows for **interactive therapies** for the whole body and comes with various equipment to make therapy more fun and engaging.



PABLO®


Lower Extremity

PABLO® Lower Extremity is a **gait analysis** and training system. It precisely assesses main gait parameters which provide essential information for treatment selection for complex walking problems.



TYMO®

TYMO® is a compact **assessment and therapy system** for the whole body. It allows to target **postural control** and **balance** in meaningful and playful manner. TYMO® comes with a variety of balance accessories for static and dynamic application in any therapy situation, e.g., standing, sitting, or supporting.



TYROSTATION

The Tyrostation is home to all individual components of PABLO® and TYMO® and provides perfect **ergonomic adaptability** for every patient.



TYROMOTION

Tyromotion Inc – USA
1341 Capital Circle, SE
Suite J
Marietta, GA 30067, USA
+ 1 678 903 8915
+ 1 888 845 4494
usa@tyromotion.com

2022 www.tyromotion.com

For more
information visit
our website or
contact us